

MONTHLY THERMOMIXER TEMPERATURE CHECK FORM

DATE _____ INDIVIDUAL(S) TESTING _____

A NIST traceable thermometer is used to check the actual temperature of the ThermoMixers.

ThermoMixer A		ThermoMixer B		ThermoMixer C		ThermoMixer D*		ThermoMixer E	
Display	Actual	Display	Actual	Display	Actual	Display	Actual	Display	Actual
56 °C		56 °C		56 °C		70 °C		56 °C	
Acceptable range 53° - 59° C (* 67° - 73° C for "ThermoMixer D" Only)									

ThermoMixer F		ThermoMixer G		ThermoMixer H		ThermoMixer I		ThermoMixer J	
Display	Actual	Display	Actual	Display	Actual	Display	Actual	Display	Actual
56 °C		56 °C		56 °C		56 °C		56 °C	
Acceptable range 53° - 59° C									

Document ID	Revision	Approval	Date Published
1801	6	Supervising Criminalist - Biology	10/23/2017 8:43:43 AM